

# SAM BalanceLab 3

- SETUP GUIDE -

for

- LAUNCH MONITOR INTEGRATION-

Trackman 4





## TrackMan Settings

- This manual refers on settings for TM4 -

### 1. Start the TrackMan software and open the settings page:

You should see a page similar to the following. Please be sure that the "Enable SDK output" checkbox is enabled:

Settings				
General	When Changing Player	Application Layout		
Players	What data do you wish to save for each player?	High Contrast Mode		
Share	All Starred None TrackMan Data:	Enable Docking in Analyze (requires restart) Reset Layouts		
About	Recorded Videos:	Power Options for Radar		
	Prompt when changing session	On battery		
	Default Ball	Dwell time after each recording 0 sec. 💙		
	Premium v	Enter power save mode after 2 min. v		
	Units	Turn off radar when leaving Live tab		
	Units Unit system (Distance, Velocity) yds/ft, mph v Normalization (convert to)	Other		
	Normalization (convert to)	✓ Save TMD files		
	Pall Promium	Delete TMD filesafter 7 days		
		Enable SDK output		
	Altitude 0 feet	Enable indoor ball flight simulation		
	Temperature 77 ° Celsius • Fahrenheit	Video Reset		
	Data	Record before impact Record after impact		
	Backup all shots, players and videos to file(s)	2,1 s 🗸 A 1,0 s 🖌 A		
	Import Tags from text file (one tag per line)	Display before impact Display after impact		
	Shot Statistic	1,5 s v ^ 0,8 s v ^		



2. Start the SAM BalanceLab 3 software and open the "Settings dialog". Open the "Launch Monitor" tab and click on Trackman.

Launch Monitor			
<ul> <li>no Launch Monitor Integration</li> <li>Flightscope</li> <li>Trackman</li> <li>Foresight API</li> <li>Foresight CSV file</li> </ul>			
IP address of Trackman computer/system: 172 . 30 . 20 . 1 Trackman running on same computer Port number: 52950			
Display units:  metric  imperial Seconds to wait for Launch Monitor impact event: 5			

#### TM4 standard IP Addresses

USB connected: address 169.254.0.1 WiFi connected: address 172.30.20.1

Enter the IP address according to your connection type and leave the port number as 52950.

Leave the checkbox "Trackman running on same computer" unchecked!

Select if you want to use metric or imperial values.

Press OK button to finish the setup.



3. Close the settings and quit the BalanceLab software. Shut down and restart both systems.

After restart of the software start a measurement in the BalanceLab software.

In the pressure curve data window select the "Launch Monitor" view.

▼ Force Plot View			
Force Plot View		-	
Launch Monitor			
975.00	4		

In the Launch Monitor View the Status should be green and say: connected to Trackman:

Launch Monitor								
Page 1 Page 2								
Carry distance	Total distance	Peak height	Club Speed	Ball Speed	Smash factor			
n/a	n/a	n/a	n/a	n/a	n/a			
Spin (total)	Spin Axis	Backspin	Sidespin	Launch vert	Classification			
n/a	n/a	n/a	n/a	n/a	n/a			
Status: - connected to FLIGHTSCOPE -								

If you now do a swing that is detected by the launch monitor system the data will be automatically be transferred to the BalanceLab software and the result numbers are displayed in the table view.

If the swing is also recorded and saved by BalanceLab the data will be saved to this data record and can be recalled in Replay and Report.



#### If the connection does not work try the following:

To enable third party applications to connect to TrackMan the SDK must be enabled in the TrackMan settings.

Enable SDK on the radar

.

Connect the radar to PC via WiFi Open a web browser In the address field, go to 172.30.20.1 Click LOGIN (no password needed) Go to SETUP and then SYSTEM Click/enable the checkbox "Legacy SDK" Click Logout

Restart Trackman system and try again to do a measurement.